

## Introduction

---

Caspase-3 is a member of caspases that plays a key role in mediating apoptosis, or programmed cell death. Upon activation, it cleaves a variety of cellular proteins, causing morphological and functional changes to cells undergoing apoptosis. The ScienCell™ Caspase-3 Assay provides a quick and convenient method to measure caspase-3 activity. The colorimetric assay is based on the spectrophotometric detection of the chromophore p-nitroanilide (pNA) after its cleavage by caspases-3 from the labeled substrate acetyl-Asp-glu-Val-Asp p-nitroanilide (Ac-DEVD-pNA). The concentration of pNA is measured by absorbance at 405 nm. The caspase-3 activity can be calculated as  $\mu\text{mol}$  of pNA released per minute per milliliter of cell lysate.

## Kit Components

| Cat. No. | # of vials | Reagent                    | Amount           | Storage     |
|----------|------------|----------------------------|------------------|-------------|
| 8228a    | 1          | Lysis Buffer               | 10 ml            | 2-8°C       |
| 8228b    | 1          | DTT Stock (1 M)            | 0.2 ml           | -20°C       |
| 8228c    | 1          | 10× Assay Buffer           | 3 ml             | 2-8°C       |
| 8228d    | 1          | Caspase-3 Substrate (2 mM) | 1 ml             | -20°C, dark |
| 8228e    | 1          | pNA Standard (40 mM)       | 20 $\mu\text{l}$ | -20°C       |

## Quality Control

---

ScienCell™ Caspase-3 Assay is applied to various concentration of active recombinant human caspase-3 (1.5 units, 5 units and 15 units per 100  $\mu\text{l}$  reaction) with ( $2 \times 10^{-3}$   $\mu\text{mol}$ ) and without caspase-3 inhibitor, according to Table 1. After incubation at 37°C for 2 hours, absorbance at 405 nm is read. Results show that the  $\text{OD}_{405\text{nm}}$  increases as the concentration of caspase-3 increases, while  $2 \times 10^{-3}$   $\mu\text{mol}$  of caspase-3 inhibitor inhibits the activity of up 15 units of caspase-3 effectively, which provides a negative control for the assay (Figure 1).

## Procedures

---

### A. Preparation of reagent

1. Aliquot and store DTT Stock (1 M) at -20°C.
2. Add appropriate volume of 1 M DTT into the Lysis Buffer to a final concentration of 5 mM (dilute 200×) before each use. Lysis Buffer with 5 mM DTT is stable for less than a week at 2-8°C.
3. Add appropriate volume of 1 M DTT into the 10× Assay Buffer to a final concentration of 50 mM (dilute 20×) before each use. 10× Assay Buffer with 50 mM DTT is stable for less than a week at 2-8°C.

### B. Preparation of pNA standard

1. Dilute 0.3 ml of 10× Assay Buffer with 50 mM DTT 10 times with DI H<sub>2</sub>O to make 3 ml of working Assay Buffer.

- Add 3  $\mu\text{l}$  of 40 mM pNA stock to 297  $\mu\text{l}$  of working Assay Buffer to make a 300  $\mu\text{l}$  solution of 400  $\mu\text{M}$  pNA.
- Obtain 8 test tubes, add 300  $\mu\text{l}$  of working Assay Buffer into each tube and label them #1 through #8.
- Add 300  $\mu\text{l}$  of the 400  $\mu\text{M}$  pNA into tube #1 and mix well to get the 200  $\mu\text{M}$  pNA standard.
- Transfer 300  $\mu\text{l}$  of 200  $\mu\text{M}$  pNA standard from tube #1 to tube #2 and mix well to get the 100  $\mu\text{M}$  pNA standard.
- Repeat step 5 for tubes #3-7 to serially dilute the pNA standards. Do not add any pNA to tube #8, which serves as the blank.
- Obtain a 96-well plate, add 100  $\mu\text{l}$ /well of each pNA standard into the 96-well plate in triplicate to generate 0.02  $\mu\text{mol}$  to  $3.125 \times 10^{-4}$   $\mu\text{mol}$ /well standard, according to the following plate format:

|   | #1                   | #2                   | #3                                 | #4                                   | #5                                    | #6                                    | #7                                     | #8    |
|---|----------------------|----------------------|------------------------------------|--------------------------------------|---------------------------------------|---------------------------------------|--|-------|
| A | 0.02 $\mu\text{mol}$ | 0.01 $\mu\text{mol}$ | $5 \times 10^{-3}$ $\mu\text{mol}$ | $2.5 \times 10^{-3}$ $\mu\text{mol}$ | $1.25 \times 10^{-3}$ $\mu\text{mol}$ | $6.25 \times 10^{-4}$ $\mu\text{mol}$ | $3.125 \times 10^{-4}$ $\mu\text{mol}$ | Blank |
| B | 0.02 $\mu\text{mol}$ | 0.01 $\mu\text{mol}$ | $5 \times 10^{-3}$ $\mu\text{mol}$ | $2.5 \times 10^{-3}$ $\mu\text{mol}$ | $1.25 \times 10^{-3}$ $\mu\text{mol}$ | $6.25 \times 10^{-4}$ $\mu\text{mol}$ | $3.125 \times 10^{-4}$ $\mu\text{mol}$ | Blank |
| C | 0.02 $\mu\text{mol}$ | 0.01 $\mu\text{mol}$ | $5 \times 10^{-3}$ $\mu\text{mol}$ | $2.5 \times 10^{-3}$ $\mu\text{mol}$ | $1.25 \times 10^{-3}$ $\mu\text{mol}$ | $6.25 \times 10^{-4}$ $\mu\text{mol}$ | $3.125 \times 10^{-4}$ $\mu\text{mol}$ | Blank |

- Read samples at 405 nm on a microtiter plate reader. Plot the standard curve of  $\text{OD}_{405\text{nm}}$  vs.  $\mu\text{mol}$  of pNA (e.g. Figure 2). Determine the equation and  $R^2$  value of the trend line.

### C. Preparation of cell lysate

- Induce apoptosis in cells by desired method.
- Harvest cell pellet for each sample. Wash the cell pellet once with PBS. Count the number of cells.
- Resuspend cells in pre chilled Lysis Buffer with 5 mM DTT at  $1 \times 10^7$  cells/100 $\mu\text{l}$ ; leave the cells on ice for 15 minutes with gentle agitation.
- Centrifuge the lysed cells at  $14,000 \times g$  in pre-cooled centrifuge for 3 minutes, transfer the supernatant to a fresh tube and discard the pellet. Cell lysate can be stored at  $-70^\circ\text{C}$  or used immediately for caspase-3 measurement.

### D. Assay procedure

- Sequentially add 20  $\mu\text{l}$  of cell lysate, 10  $\mu\text{l}$  of 10 $\times$  Assay Buffer with 50 mM DTT, 60  $\mu\text{l}$  of DI  $\text{H}_2\text{O}$  and 10  $\mu\text{l}$  of 2 mM Caspase-3 Substrate to each well of a 96 well plate. Prepare a couple of blank wells by mixing 20  $\mu\text{l}$  of Lysis Buffer with 5 mM DTT, 10  $\mu\text{l}$  of 10 $\times$  Assay Buffer with 50 mM DTT, 60  $\mu\text{l}$  of DI  $\text{H}_2\text{O}$  and 10  $\mu\text{l}$  of 2 mM Caspase-3 Substrate in each well of the 96 well plate. Incubate at  $37^\circ\text{C}$  for 2-4 hours or until a yellowish color is developed. Record the time of reaction in minutes.
- Read samples at 405 nm on a microtiter plate reader.

### E. Calculation

- Subtract the averaged  $\text{OD}_{405\text{nm}}$  of the blank wells from each of the sample well to get the calibrated  $\text{OD}_{405\text{nm}}$  values of the sample wells. Suppose the equation of the trend line of the pNA standard

curve is  $y = Ax + B$ , calculate the  $\mu\text{mol}$  of the pNA released in each sample well as follows:

$$pNA = \frac{OD_{405nm} - B}{A}$$

2. Calculate the caspase-3 activity in  $\mu\text{mol}$  pNA released per min per ml of cell lysate as follows:

$$Activity, \mu\text{mol pNA /min /mL} = \frac{\mu\text{mol of pNA}}{0.1\text{mL (lysate volume)} \times t}$$

Where t is the reaction time in minutes.

Table 1. Reaction scheme for caspase-3 positive control with and without inhibitor.

|  | <b>Caspase-3<br/>1 unit/<math>\mu</math>l</b> | <b>Caspase-3 Substrate<br/>2 mM</b> | <b>Caspase-3 Inhibitor<br/>0.2 mM</b> | <b>10<math>\times</math> Assay Buffer<br/>with 50 mM DTT</b> | <b>DI H<sub>2</sub>O</b> |
|--|---|-------------------------------------|---------------------------------------|--|--------------------------|
| Caspase-3<br>positive control<br>without inhibitor | 1.5 $\mu$ l                                   | 10 $\mu$ l                          | --                                    | 10 $\mu$ l   | 78.5 $\mu$ l             |
|  | 5 $\mu$ l                                     | 10 $\mu$ l                          | --                                    | 10 $\mu$ l   | 75 $\mu$ l               |
|  | 15 $\mu$ l                                    | 10 $\mu$ l                          | --                                    | 10 $\mu$ l   | 65 $\mu$ l               |
| Caspase-3<br>positive control<br>with inhibitor    | 1.5 $\mu$ l                                   | 10 $\mu$ l                          | 10 $\mu$ l                            | 10 $\mu$ l   | 68.5 $\mu$ l             |
|  | 5 $\mu$ l                                     | 10 $\mu$ l                          | 10 $\mu$ l                            | 10 $\mu$ l   | 65 $\mu$ l               |
|  | 15 $\mu$ l                                    | 10 $\mu$ l                          | 10 $\mu$ l                            | 10 $\mu$ l   | 55 $\mu$ l               |




Figure 1. ScienCell™ Caspase-3 Assay kit applied to various amount of caspase-3 positive control with and without inhibitor (Table 1).




Figure 2. Standard curve of OD<sub>405nm</sub> vs. pNA in μmol.